



# CROWN REMOVAL

IN **2** MINUTES



## CLINICAL CASES

AVAILABLE ON OUR WEBSITE

[www.wamkey.com](http://www.wamkey.com)



Bridge removal



Drilling window and thin tunnel above each prep.



Controlling the depth of each tunnel (must reach the center)



Inserting and turning the working tip in each tunnel breaks the seal.



Removing the crown gives a better access for RCT



The bridge can be reused at least as a temporary



# WAMkey

## EFFICIENCY

TESTED & SUPERIOR TO

# 90 %\*

- Fast & simple**
- Atraumatic** for the supporting tooth
- Saves drills and preserves** handpieces
- The crown (or the bridge) can be **re-used**

*\* Evaluation Conducted by Dental Advisor in 2012, on 320 crowns (all types of crowns and sealings)*



**WAM**

Tech'Indus B21  
645 rue Mayor de Montricher  
13854 Aix-en-Provence Cedex 3

TEL

+33 (0)4 42 26 80 64

+33 (0)4 13 41 57 87

    
[www.wamkey.com](http://www.wamkey.com)